

Freestyle Flex™ Double Electric Connected Breast Pump

101037980



The Freestyle Flex™ Breast Pump is a double electric daily-use pump designed for the on-the-go mom who pumps several times a day. Designed for her active lifestyle, the Freestyle Flex is light, compact, and has a USB rechargeable lithium-ion battery to give mom the flexibility to express wherever and whenever she likes – with no compromise on pumping performance. Weighing less than 1 pound, mom can simply slip it in her purse or baby bag and she's ready to get on with her busy day.

Breast Pump Features

Reliable Closed System with Overflow Protection
prevents breast milk from entering tubing

Compact size & rechargeable battery
provides the flexibility to express wherever and whenever mom chooses

Bluetooth® connectivity to MyMedela® app
helps mom easily track and record pumping sessions

Warranty
Warranted to the original participant to be free from defects in material and workmanship for the period of one year for pump mechanism (90 days for parts and accessories).

Regulatory Info
All components that come in contact with mother's milk are made without BPA and meet US Consumer Product Safety Improvement Act (CPSIA 2008), CA AB 1108 requirements for DEHP and other phthalates for all components that contact mother's milk.



The combination of double pumping with 2-Phase Expression® Technology is proven to produce more milk in less time, when compared to single pumping**

Get connected



Freestyle Flex connects to the **MyMedela® app**, to give mom additional support and encouragement. She can use the app to track and log her pumping sessions, and create a personalized profile tailored to her, her baby, and feeding choices.

Intuitive multifunctional interface with:

- Count-up timer
- Delay start & pause functions
- Pump setting & battery charge indicators

USB cord for powering and charging anywhere you are*

Detachable lanyard loops around mom's wrist to keep breast pump secure

Built-in rechargeable battery for up to 2 hours of pumping on a full charge - no need to search for an outlet.

Compact design fits easily into mom's palm, pocket, or bag, and weighs less than 1 pound

*Only use the power adaptor and cable provided with the breast pump.

Retail Product Specifications:

Item #	Description	UPC Code	MSRP	Unit Count	Unit Dimensions (H x W x D) Inches	Unit Weight (lbs)	Qty (per case)	Shipping Carton Dimensions (HxWxD) inches	Carton Weight (lbs)
101037980	Freestyle Flex™ Breast Pump	0 20451 37980 7	\$379.99	1	9.45 x 13.43 x 7.05	5.09	1	10.43 x 15.55 x 7.76	5.86
101038442	Freestyle Flex™ Power Adaptor	0 20451 38442 9	\$29.99	1	3.35 x 3.35 x 2.64	0.18	4	2.95 x 7.09 x 7.09	1.00
101038234	Freestyle Flex™ Tubing	0 20451 38234 0	\$19.99	1	6.20 x 4.25 x 1.50	0.12	6	5.83 x 9.60 x 4.72	1.00
101038233	PersonalFit Flex™ Connectors	0 20451 38233 3	\$29.99	2	5.44 x 5.28 x 2.40	0.31	6	8.90 x 11.02 x 7.64	2.20

To order or for more information, contact Medela Customer Service at customer.service@medela.com or 888-MEDELA8 (888-633-3528).

Medela LLC, 1101 Corporate Drive, McHenry, IL 60050
 Phone: (800) 435-8316 or (815) 363-1166 Email: customer.service@medela.com www.medela.com

Content appearance may vary from pictures.

References: *Prime DK et al. Breastfeed Med. 2012; 7(2):100-106. *Prime DK et al. 6th ABM Europe Conference, Rotterdam, NL; 2018.

Medela word mark and logo, 2-Phase Expression and MyMedela are registered in the U.S. Patent and Trademark Office. Freestyle Flex is a trademark of Medela. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks is used under license.

1549493 A 1119 © 2019 Medela 4551309